

## D. Draping the Point Setter

Drape **before attaching** the holder or endoscope to the Point Setter.

### 1. Drapes for the Point Setter



**Drapes CANNOT BE REUSED**


Drapes are disposable. Keep adequate stock.  
To Order drapes see Mitaka USA below (Catalog #4667) or your local Karl Storz distributor.

### 2. Open the bag



Open the bag from the side as indicated by the arrow.  
Be careful that good sterile technique is used at all times.

### 3. Cover the Point Setter with the drape.



Secure the drape with adhesive tape, which is attached to the drape. When securing the drape, do not tighten it excessively.

## E. Attaching the Endoscope to the Point Setter

### 1. Installing an endoscope



After draping the Point Setter, attach the endoscope holder to the Point Setter. Secure firmly so no movement is present.

Now the endoscope may be attached to the endoscope holder.

Use the appropriate endoscope holder for each different type of endoscope.

**NOTE:** When installing the endoscope holder over the draped Point Setter, be careful not to tear the drape.

**NOTE:** The surgeon may choose to connect the endoscope before or during surgery.

### 2. Adjustment of manual ball joint



For optimal positioning, adjust the angle with the manual ball joint adjust knob.

### 3. Safety lock knob



**CAUTION!**

During operation, the **safety lock knob** must remain in the "LOCK" position except when moving the arm.

**NOTE:** When the safety lock knob is in the "Lock" position, it prevents accidental release of the arm.

#### 4. Ready to use



Pushing the free switch (which is colored gold) allows the arm to move freely.

Letting go of free switch should lock brakes after a very short time.



#### CAUTION!

The holding force of the Point Setter depends on the position of the device which you are going to install. The farther away from connection the less weight you can hold (due to leverage).

Make sure the device has enough holding force before operation.

The further an accessory is attached from the Point Setter, the smaller the weight that can be safely held. Careful following of this knowledge is critical for safe use of Point Setter.

